



MRI Information Insert

Non-clinical testing has demonstrated the device is MR Conditional. It can be scanned safely under:

- static magnetic field of 3 Tesla or less
- spatial gradient field of 330 Gauss/cm or less
- maximum specific absorption rate (SAR) of 6 W/kg for 30 minutes of scanning.

In non-clinical testing, the device produced a temperature rise of less than 0.5 °C at a maximum specific absorption rate (SAR) of 6 W/kg for 30 minutes of MR scanning in a 3T Siemens Trio with software version VA25.

For Minimal Image artifact

- MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

* Bard, PowerPort and the color purple are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

© Copyright 2007 C. R. Bard, Inc. All rights reserved.

Bard Access Systems, Inc.

Salt Lake City, UT 84116 USA 801-595-0700
Clinical Information Hotline: 800-443-3385
Ordering Information: 800-545-0890
www.bardaccess.com
www.portadvantage.com
www.powerportadvantage.com

BARID