

Aspira^{*}
DRAINAGE SYSTEM

Compassionate Treatment



BARD
ACCESS SYSTEMS

*Advancing the Delivery of Health Care.**

Do Your Patients Require *Frequent Thoracentesis or Paracentesis?*

Bard now offers a compassionate home drainage option that is easy to place and simple to use.

) **Quality** (

High quality kit components, including catheter with Hydro-Glide* coated stylet to decrease leakage during insertion.


) **Simplicity** (

The large 1L volume drainage collection bag is simple for a patient or caregiver to use.


) **Options** (

Multiple packaging options provide the ideal solutions for distribution of patient drainage supplies.


) **PLACEMENT OPTIONS** (



Antegrade



Retrograde



Over-the-Wire

Aspira*

DRAINAGE SYSTEM

Low Vacuum Siphon
Activating Pump


1L Drainage Bag

Attachable Valve

15.5 Fr. Silicone
Catheter

Cuff

Catheter with
Hydro-Glide* Coated Stylet



) ORDERING INFORMATION (

Everything you need in **ONE** package.

Aspira Pleural Drainage System*

Product Code	Product Code Description
4991507	Aspira* Pleural Drainage Complete Kit Each Kit Includes: 1 Insertion Tray with Pleural Catheter 5 Pleural Drainage Kits 5 Dressing Kits Voucher for 20 Drainage Kits and 20 Dressing Kits
4991508	Aspira* Pleural Drainage Essential Kit Each Kit Includes: 1 Insertion Tray with Pleural Catheter 20 Drainage Kits 20 Dressing Kits

Aspira Drainage Accessories*

Product Code	Product Code Description
4991502	1,000 ml Aspira* Drainage Bag Each Kit Includes: 5 Drainage Kits
4991503	Dressing Kit Each Kit Includes: 5 Dressing Kits
4991505	Luer Adapters 5 Ea./Box
4991506	Repair Kit 5 Ea./Box

Aspira Peritoneal Drainage System*

Product Code	Product Code Description
4991207	Aspira* Peritoneal Drainage Complete Kit Each Kit Includes: 1 Insertion Tray with Peritoneal catheter 5 Drainage Kits 5 Dressing Kits Voucher for 20 Drainage Kits and 20 Dressing Kits
4991208	Aspira* Peritoneal Drainage Essential Kit Each Kit Includes: 1 Insertion Tray with Peritoneal catheter 20 Drainage Kits 20 Dressing Kits

Aspira Patient Education*

Product Code	Product Code Description
4991555	Pleural Patient Guide DVD
4991556	Pleural Patient Guide VHS
	Peritoneal Video (Available Soon)

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Learn More - For more information visit: www.myaspira.com



ACCESS SYSTEMS

Bard Access Systems, Inc.

605 North 5600 West
 Salt Lake City, Utah 84116 USA
 801-522-5000
www.bardaccess.com

Customer Service: 800-545-0890
 Clinical Information: 866-893-2691
 Home Supply Ordering: 866-443-3471