

Full Nursing Tray components:

- 1 PowerPICC catheter with Pre-loaded stylet
- 2 Anchor, Tape Strips
- 1 ChloraPrep[®] One-Step Applicator
- 1 Drape, Absorbent
- 1 Drape, Fenestrated
- 1 Filter Straw[®], 18 Ga. x 1.5 in. (per lumen)
- 1 Gloves
- 1 Lidocaine HCl 1%, 5 cc Ampule
- 2 Mask without Shield
- 1 SafetyGlide Needle, 25 Ga. x 5/8 in.
- 1 Scissors, Disposable
- 1 Sodium Chloride (Saline) 0.9%, 10 cc Ampule / Lumen
- 1 Sorbaview[®] Dressing, Transparent
- 6 Gauze Pads: 4 in. x 4 in.
- 10 Gauze Pads: 2 in. x 2 in.
- 1 Statlock[®] Securement Device
- 1 Surgical Tape
- 3 Syringes, Luer Lock
- 1 Tape Measure, Paper
- 1 Tourniquet
- 1 Towel, Absorbent
- 1 Injection Cap (per lumen)
- 1 Wipe, Alcohol

Microintroducer contents :

- 1 Scalpel, Safety #11
- 1 Guidewire, 0.018 in x 35 cm
- 1 Introducer, IV Catheter, Safety, 20 Ga. x 1.75 in.
- 1 Introducer, Safety Seldinger Needle 21 Ga. x 2.75 in.
- 1 MicroEZ Microintroducer

Basic Nursing Tray Components

- 1 PowerPICC catheter with Pre-loaded stylet
- 1 Injection Cap (per lumen)
- 1 Connector, T-Lock with Injection Port
- 1 StatLock[®] Securement Device

Microintroducer contents:

- 1 Scalpel, Safety #11
- 1 Guidewire, 0.018 in x 35 cm
- 1 Introducer, IV Catheter, Safety, 20 Ga. x 1.75 in.
- 1 Introducer, Safety Seldinger Needle 21 Ga. x 2.75 in.
- 1 MicroEZ Microintroducer

Nursing	Configuration	MicroEz MicroIntroducer Length	Nursing Guidewire Length
	5 Fr SL	5 cm	35 cm
	5 Fr DL	5 cm	35 cm
	6 Fr DL	10 cm	45 cm
	6 Fr TL	10 cm	35 cm

Interventional Radiology Tray Components

- 1 PowerPICC Catheter
- 1 Gold Tungsten Tip Guidewire, Marked, .018"
- 1 PTFE Peel-Apart Introducer and Dilator
- 1 Injection Cap (per lumen)
- 1 Introduction Needle, 21ga x 7cm length
- 1 StatLock® Securement Device
- 1 Scalpel, #11
- 1 Syringe, 12cc
- 1 Tape Measure

IR	Configuration	Microintroducer Length	Interventional Radiology Guidewire Length	Guidewire configuration options
	5 Fr SL	10 cm	70 cm, 135 cm	Standard, Reinforced Tip, Nitinol
	5 Fr DL	10 cm	70 cm, 135 cm	Standard, Reinforced Tip, Nitinol
	6 Fr DL	10 cm	70 cm, 135 cm	Standard, Reinforced Tip, Nitinol
	6 Fr TL	10 cm	70 cm, 135 cm	Standard, Reinforced Tip, Nitinol