


PowerPICC* Triple Lumen Catheter: The ULTIMATE CHOICE for Access


Now Available...

5 Fr. Triple-Lumen PowerPICC* Catheter
A Smaller PowerPICC* Catheter Option



BAIRD
ACCESS SYSTEMS

*Advancing the Delivery of Health Care.**



The best-selling PowerPICC* French size with the added utility of a third lumen to provide peripheral vascular access to a wider range of patients.

Benefits Include:

- 5 ml/sec max flow rate
- Indicated for CVP monitoring
- Designed for multiple infusion therapies as well as blood sampling
- Center Lumens designed to accommodate high flow rates (999 ml/hr pump flow rate)
- Long term peripheral access
- Provides vascular access for patients in the ICU and after transfer out of the ICU
- Lower placement risks (compared to Acute CVCs) - due to peripheral insertion in upper arm
- Can be placed at the bedside by a skilled nurse

Ordering Information

Product Code	Product Description	Tray	Lumen	Gravity Flow	Priming Volume	Power Injection Flow Rate	
5 Fr.	3385108	5 Fr. PowerPICC* Triple Lumen catheter	Full Nursing Tray	18 Ga 20 Ga	982 ml/hr 131 ml/hr	0.76 ml 0.43 ml	5 ml/sec N/A
	9385108	5 Fr. PowerPICC* Triple Lumen catheter w/ Sherlock* Stylet	Full Nursing Tray	18Ga 20Ga	982 ml/hr 131 ml/hr	0.76 ml 0.43 ml	5 ml/sec N/A
	3385108D	5 Fr. PowerPICC* Triple Lumen catheter w/ Max Barrier	Full Max Barrier Nursing Tray	18 Ga 20 Ga	982 ml/hr 131 ml/hr	0.76 ml 0.43 ml	5 ml/sec N/A
	9385108D	5 Fr. PowerPICC* Triple Lumen catheter w/ Max Barrier & Sherlock* Stylet	Full Max Barrier Nursing Tray	18 Ga 20 Ga	982 ml/hr 131 ml/hr	0.76 ml 0.43 ml	5 ml/sec N/A
	3385118	5 Fr. PowerPICC* Triple Lumen catheter	Basic Nursing Tray	18 Ga 20 Ga	982 ml/hr 131 ml/hr	0.76 ml 0.43 ml	5 ml/sec N/A
	9385118	5 Fr. PowerPICC* Triple Lumen catheter w/ Sherlock* Stylet	Basic Nursing Tray	18 Ga 20 Ga	982 ml/hr 131 ml/hr	0.76 ml 0.43 ml	5 ml/sec N/A
	3385335	5 Fr. PowerPICC* Triple Lumen catheter - 135cm Nitinol Guidewire	Basic IR OTW Tray	18 Ga 20 Ga	982 ml/hr 131 ml/hr	0.76 ml 0.43 ml	5 ml/sec N/A
	3385355	5 Fr. PowerPICC* Triple Lumen catheter - 70cm Nitinol Guidewire	Basic IR Tray	18 Ga 20 Ga	982 ml/hr 131 ml/hr	0.76 ml 0.43 ml	5 ml/sec N/A
6 Fr.	3386108	6 Fr. PowerPICC* Triple Lumen catheter	Full Nursing Tray	17 Ga 19 Ga	1163 ml/hr 275 ml/hr	0.76 ml 0.47 ml	5 ml/sec N/A
	9386108	6 Fr. PowerPICC* Triple Lumen catheter w/ Sherlock* Stylet	Full Nursing Tray	17 Ga 19 Ga	1163 ml/hr 275 ml/hr	0.76 ml 0.47 ml	5 ml/sec N/A
	3386108D	6 Fr. PowerPICC* Triple Lumen catheter w/ Max Barrier	Full Max Barrier Nursing Tray	17 Ga 19 Ga	1163 ml/hr 275 ml/hr	0.76 ml 0.47 ml	5 ml/sec N/A
	9386108D	6 Fr. PowerPICC* Triple Lumen catheter w/ Max Barrier & Sherlock* Stylet	Full Max Barrier Nursing Tray	17 Ga 19 Ga	1163 ml/hr 275 ml/hr	0.76 ml 0.47 ml	5 ml/sec N/A
	3386105	6 Fr. PowerPICC* Triple Lumen catheter	Basic Nursing Tray	17 Ga 19 Ga	1163 ml/hr 275 ml/hr	0.76 ml 0.47 ml	5 ml/sec N/A
	9386105	6 Fr. PowerPICC* Triple Lumen catheter w/ Sherlock* Stylet	Basic Nursing Tray	17 Ga 19 Ga	1163 ml/hr 275 ml/hr	0.76 ml 0.47 ml	5 ml/sec N/A
	3386335	6 Fr. PowerPICC* Triple Lumen catheter - 135cm Nitinol Guidewire	Basic IR OTW Tray	17 Ga 19 Ga	1163 ml/hr 275 ml/hr	0.76 ml 0.47 ml	5 ml/sec N/A
	3386355	6 Fr. PowerPICC* Triple Lumen catheter - 70cm Nitinol Guidewire	Basic IR Tray	17 Ga 19 Ga	1163 ml/hr 275 ml/hr	0.76 ml 0.47 ml	5 ml/sec N/A

The maximum recommended infusion rate is 5 ml/sec for power injection of contrast media. The maximum pressure of power injectors used with the PowerPICC* catheter may not exceed 300psi.

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Learn More - To learn more about the PowerPICC* Triple Lumen Catheter and how it can help streamline your hospital's workflow, call Bard Access Systems at 1-800-545-0890 or visit us on the web at www.bardaccess.com.