




PowerPICC
solo Catheter

Power. Saline. One.*

BAIRD
ACCESS SYSTEMS

*Advancing the Delivery of Health Care.**



PowerPICC SOLO* catheter is a universal catheter solution designed for IV fluid and blood administration, power injection, and saline care and maintenance.

Power. Saline. One.*

Power.


The Power of Purple*

The Efficiency of Power Injection for Contrast Enhanced CT Scans and PICC Therapy.

- Maximum Injection Rates - allows injection of contrast media for CECT scans at a maximum rate of **5 ml/sec**
- Easy Identification - catheters are clearly labeled to differentiate them from all non-power injectable PICCs
- CVP Monitoring - all PowerPICC SOLO* catheters are indicated for central venous pressure monitoring
- Reverse Taper Design - gently plugs the insertion site and provides superior kink resistance
- Patient Satisfaction - improves patient comfort due to decreased needle sticks

THE POWER OF PURPLE*

Power Injectable PICC



PowerPICC
solo Catheter

Saline.

The Simplicity of Saline*



The Simplicity of Saline Care and Maintenance Through Innovative Valve Technology.

- **Simplified Care and Maintenance** - recommended weekly maintenance is reduced to a single saline flush
- **Continuum of Care** - same flushing protocol for hospital and home health
- **Fewer Complications** - limits risk associated with heparin use (HIT)
- **Reduced Blood Reflux** - compared to open-ended catheters
- **Reduced Costs** - decreases PICC associated heparin costs
- **Clamp Free Infusion Therapy** - valve eliminates the need for clamps and the maintenance challenges associated with needless injection caps
- **Superior Clearance** - valve clears after recommended saline flush



Infusion



Aspiration

BLUE COLOR

Symbolizes **Saline**
Care and Maintenance

1

One.

One Valve. One PICC.*

The Versatility of **One** PICC Designed to Meet a Variety of Clinician and Patient Needs.

- **Simplified Flushing** - recommended weekly maintenance is reduced to a single flush
- **Continuum of Care** - standardizes care and maintenance for hospital and home health
- **Broad Vascular Access Needs** - designed for multiple clinical conditions across the hospital
- **Multiple Indications** - allows IV fluid and blood infusion, as well as power injection, and central venous pressure monitoring
- **Many Configurations** - 4 and 5 Fr. Single-Lumen, 5 and 6 Fr. Dual-Lumen, 6 Fr. Triple-Lumen
 - **Variety of Placement Technologies** -
 - Nursing** - preloaded Sherlock* and non-Sherlock* stylets
 - Interventional Radiology** - allows for over-the-guidewire placement

5 Fr. Dual-Lumen
PowerPICC SOLO* Catheter
(Also Available in 6 Fr.)

6 Fr. Triple-Lumen
PowerPICC SOLO* Catheter

4 Fr. Single-Lumen
PowerPICC SOLO* Catheter
(Also Available in 5 Fr.)





Power. Saline. One.*

For more information visit: www.powerpiccsolo.com

BARD

ACCESS SYSTEMS

Bard Access Systems, Inc.

Salt Lake City, Utah 84116
801-595-0700
www.bardaccess.com

Customer Service: 800-545-0890
Clinical Information: 800-443-3385

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

*Bard, PowerPICC, PowerPICC SOLO, The Simplicity of Saline, One Valve. One PICC., Power. Saline. One., The Power of Purple, Sherlock, "Advancing the Delivery of Health Care." and the color Purple are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.
© Copyright 2007 C. R. Bard, Inc. All rights reserved.

MC-0130-00 0710R