

## New and Improved Kit Components



Power.  
The Power of Purple\*



- Kink resistant guidewire
- Bendable proximal end
- Length of guidewire offers greater placement versatility



Saline.  
The Simplicity of Saline\*

### Product Specifications:

Type of Catheter	Product Code	Tray	Gravity Flow (ml/hr)	Priming Volume (ml)	Gauge	Introducer Handle Color
4 Fr. SL PowerPICC SOLO* Catheter	3194118 3194108	Basic w/Microintroducer Full w/Microintroducer	1024	0.73	18	
5 Fr. SL PowerPICC SOLO* Catheter	3195118 3195108	Basic w/Microintroducer Full w/Microintroducer	978	0.73	18	
5 Fr. DL PowerPICC SOLO* Catheter	3295118 3295108	Basic w/Microintroducer Full w/Microintroducer	498/498	0.65/0.65	18/18	
6 Fr. DL PowerPICC SOLO* Catheter	3296118 3296108	Basic w/Microintroducer Full w/Microintroducer	740/740	0.75/0.75	18/18	
6 Fr. TL PowerPICC SOLO* Catheter	3396118 3396108	Basic w/Microintroducer Full w/Microintroducer	R:560 G:308 W:308	R:0.71 G:0.57 W:0.57	R:17 G:19 W:19	

One.  
One Valve. One PICC.\*

### Product Specifications with Sherlok\* Stylet

Type of Catheter	Product Code	Tray	Gravity Flow (ml/hr)	Priming Volume (ml)	Gauge	Introducer Handle Color
4 Fr. SL PowerPICC SOLO* Catheter	9194118 9194108	Basic w/Microintroducer Full w/Microintroducer	1024	0.73	18	
5 Fr. SL PowerPICC SOLO* Catheter	9195118 9195108	Basic w/Microintroducer Full w/Microintroducer	978	0.73	18	
5 Fr. DL PowerPICC SOLO* Catheter	9295118 9295108	Basic w/Microintroducer Full w/Microintroducer	498/498	0.65/0.65	18/18	
6 Fr. DL PowerPICC SOLO* Catheter	9296118 9296108	Basic w/Microintroducer Full w/Microintroducer	740/740	0.75/0.75	18/18	
6 Fr. TL PowerPICC SOLO* Catheter	9396118 9396108	Basic w/Microintroducer Full w/Microintroducer	R:560 G:308 W:308	R:0.71 G:0.57 W:0.57	R:17 G:19 W:19	

To learn more about how the PowerPICC SOLO\* Catheter can help your hospital, call Bard Access Systems at **800-545-0890** or visit us on the web at [www.powerpiccso.com](http://www.powerpiccso.com) or [www.bardaccess.com](http://www.bardaccess.com).



Power. Saline. One.\*

**BARD**

ACCESS SYSTEMS

**Bard Access Systems, Inc.**

Salt Lake City, Utah 84116  
801-595-0700  
[www.bardaccess.com](http://www.bardaccess.com)

Customer Service: 800-545-0890  
Clinical Information: 800-443-3385

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

\*Bard, PowerPICC, PowerPICC SOLO, The Simplicity of Saline, One Valve, One PICC., Power, Saline, One., The Power of Purple, MicroEZ, Flexura, Sherlock and the color Purple are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.  
© Copyright 2008 C. R. Bard, Inc. All rights reserved.

MC-0132-01 0804R